

stryker



Driving
performance,
**delivering
versatility**

**Malis and SILVERGlide
Non-Stick Bipolar Forceps**

Your **superior experience** starts with our people

Supporting you throughout your procedures

People

Our people stand out. They're dedicated, highly skilled experts that care deeply and listen carefully. When you need a partner in improving care, we're by your side.

Culture

We create a dynamic, supportive environment that enables amazing people to do amazing things. Count on us to do what's right and deliver—every time.

Product innovation

Innovation is in our DNA—from Dr. Homer Stryker to the researchers, engineers and surgeons creating solutions that are reshaping neurotechnology today.

Medical education

Innovation and skill come together in our educational outreach to providers. Each year, we help surgical teams around the world master the latest medical technologies and care practices.

Research and development


We amplify our investments in R&D by listening to our customers and collaborating with them to create solutions that are simply superior.

Outcomes

We share your values and your commitment to innovative care that heals effectively and efficiently. With you, we're driven to make healthcare better.



All of our non-stick bipolar forceps **resist tissue buildup** and provide **excellent visualization**. Regardless of your preference for single-use Malis or reusable SILVERGlide, you can expect **precision, control** and a responsive extension of your hands and eyes from us.



Disposable Malis

Standard Bipolar Forceps

Flexible lightweight cord and remote connection plug helps achieve balance

Spring force elevates control and tactile feedback

Single use ensures excellent tip alignment and coagulation performance

Sleek bayonet design improves visualization

Non-stick technology enhances performance



7" (18cm)

6760-180-005 0.5mm tip
6760-180-010 1.0mm tip
6760-180-015 1.5mm tip

8" (20cm)

6760-200-005 0.5mm tip
6760-200-010 1.0mm tip
6760-200-015 1.5mm tip

9" (23cm)

6760-230-005 0.5mm tip
6760-230-010 1.0mm tip
6760-230-015 1.5mm tip

Dual Irrigating Bipolar Forceps

Transparent fine tubing allows for better observation

Dual irrigation design maintains moisture of the tips and vessels



7" (18cm)

6770-180-005 0.5mm tip
6770-180-010 1.0mm tip
6770-180-015 1.5mm tip

8" (20cm)

6770-200-005 0.5mm tip
6770-200-010 1.0mm tip
6770-200-015 1.5mm tip

9" (23cm)

6770-230-005 0.5mm tip
6770-230-010 1.0mm tip
6770-230-015 1.5mm tip

Slim Bipolar Forceps

Slim profile provides access in narrow fields



8" (20cm)


6780-200-005 0.5mm tip
6780-200-010 1.0mm tip
6780-200-015 1.5mm tip

9" (23cm)

6780-230-005 0.5mm tip
6780-230-010 1.0mm tip
6780-230-015 1.5mm tip

SILVERGlide Non-Stick

ProSeries Bipolar Forceps



Standard two-pin round connector works with most standard foot-switching bipolar electro-surgical generators

Silver plating reduces collateral heating, "hot spots" and sticking for efficient performance

Ultra-thin tines allow better access in surgery

Lightweight handle for easy manipulation and balance

Guide-stop minimizes scissoring and increases control



5.7" (14.5cm) bayonet

- 6730-145-004 0.4mm
- 6730-145-007 0.7mm
- 6730-145-010 1.0mm
- 6730-145-012 1.2mm
- 6730-145-015 1.5mm

6.3" (16cm) bayonet 6730-160-004

- 0.4mm 6730-160-007 0.7mm
- 6730-160-010 1.0mm 6730-160-012 1.2mm
- 6730-160-504, 45° up 0.4mm

7.5" (19cm) bayonet

- 6730-190-004 0.4mm
- 6730-190-007 0.7mm
- 6730-190-010 1.0mm
- 6730-190-012 1.2mm
- 6730-190-015 1.5mm
- 6730-190-020 2.0mm
- 6730-190-504, 45° up 0.4mm
- 6730-190-512, 45° up 1.2mm
- 6730-190-612, 45° down 1.2mm

8.6" (22cm) bayonet

- 6730-220-004 0.4mm
- 6730-220-007 0.7mm
- 6730-220-012 1.2mm
- 6730-220-015 1.5mm
- 6730-220-020 2.0mm
- 6730-220-107, 15° up 0.7mm
- 6730-220-307, 30° up 0.7mm
- 6730-220-507, 45° up 0.7mm
- 6730-220-607, 45° down 0.7mm
- 6730-220-512, 45° up 1.2mm
- 6730-220-612, 45° down 1.2mm

7.5" (19cm) straight bayonet

- 6735-190-007 0.7mm
- 6735-190-010 1.0mm
- 6735-190-310, 30° up 1.0mm
- 6735-190-320, 30° up 2.0mm



9.4" (24cm) bayonet

- 6730-240-004 0.4mm
- 6730-240-007 0.7mm
- 6730-240-012 1.2mm
- 6730-240-015 1.5mm
- 6730-240-020 2.0mm
- 6730-240-507, 45° up 0.7mm
- 6730-240-607, 45° down 0.7mm
- 6730-240-807, 80° up 0.7mm
- 6730-240-907, 80° down 0.7mm



11" (28cm) bayonet

- 6730-280-012 1.2mm
- 6730-280-015 1.5mm

Keyhole Maximum Access Bipolar Forceps



5.5" (14cm) keyhole bayonet

- 6750-140-004 0.4mm
- 6750-140-007 0.7mm

6.7" (17cm) keyhole bayonet

- 6750-170-004 0.4mm
- 6750-170-007 0.7mm

7.8" (20cm) keyhole bayonet

- 6750-200-004 0.4mm
- 6750-200-007 0.7mm
- 6750-200-012 1.2mm
- 6750-200-015 1.5mm

9.2" (23cm) keyhole bayonet

- 6750-230-007 0.7mm
- 6750-230-012 1.2mm
- 6750-230-015 1.5mm
- 6750-230-107, 15° up 0.7mm
- 6750-230-307, 30° up 0.7mm

10.2" (25cm) keyhole bayonet

- 6750-250-007 0.7mm
- 6750-250-012 1.2mm
- 6750-250-015 1.5mm
- 6750-250-507, 45° up 0.7mm
- 6750-250-607, 45° down 0.7mm
- 6750-250-807, 80° up 0.7mm
- 6750-250-907, 80° down 0.7mm

11.7" (28cm) keyhole bayonet

- 6750-280-012 1.2mm
- 6750-280-015 1.5mm

SILVERglide Non-Stick

Precision Bipolar Forceps



6.1" (15.5cm) SP bayonet

6720-155-004 0.4mm
6720-155-007 0.7mm
6720-155-012 1.2mm
6720-155-015 1.5mm

8.7" (22cm) SP bayonet

6720-220-004 0.4mm
6720-220-007 0.7mm
6720-220-012 1.2mm
6720-220-015 1.5mm
6720-220-020 2.0mm

7.0" (18cm) SP straight

6725-180-007 0.7mm
6725-180-010 1.0mm
6725-180-410, 30° up 1.0mm
6725-180-420, 30° up 2.0mm

9.7" (24cm) SP bayonet

6720-240-004 0.4mm
6720-240-007 0.7mm
6720-240-015 1.5mm
6720-240-020 2.0mm

7.5" (19cm) SP bayonet

6720-190-004 0.4mm
6720-190-007 0.7mm
6720-190-010 1.0mm
6720-190-015 1.5mm
6720-190-020 2.0mm
6720-190-512, 45° up 1.2mm
6720-190-612, 45° down 1.2mm

Standard Bipolar Forceps



4.0" (10cm) jewelers

6715-100-004 0.4mm
6715-100-007 0.7mm

7.5" (19cm) cushioning

6710-190-004 0.4mm
6710-190-007 0.7mm
6710-190-010 1.0mm
6710-190-012 1.2mm
6710-190-015 1.5mm
6710-190-020 2.0mm

5.0" (12.5cm) adson

6715-125-004 0.4mm
6715-125-007 0.7mm
6715-125-012 1.2mm

7.8" (20cm) scoville-greenwood

6710-200-015 1.5mm

6.0" (15cm) semkin

6715-150-012 1.2mm

8.3" (21cm) cushioning

6710-210-007 0.7mm
6710-210-012 1.2mm
6710-210-015 1.5mm
6710-210-020 2.0mm

Irrigating Bipolar Forceps

Hollow tine design
brings irrigation to the
tip of the forceps



7.5" (19cm) bayonet

6740-190-007 0.7mm
6740-190-012 1.2mm

8.7" (22cm) bayonet

6740-220-007 0.7mm
6740-220-012 1.2mm





10.2" (25cm) bayonet

6740-250-007 0.7mm
6740-250-012 1.2mm

Separate fluid and
electrical connectors
accept a variety of cord-tubing sets

What sets SILVERGlide apart?

The elimination of tissue buildup allows you to work without interruption. These non-stick characteristics continue to perform use after use. The SILVERGlide alloy's level of thermal diffusivity is 37 times higher than that of stainless steel commonly used in forceps.

	Standard forceps	SILVERGlide Non-Stick	
"Hot spots" from sparks occur during coagulation			"Hot spots" start to occur during coagulation
The "hot spots" cause localized sticking			Heat is diffused rapidly to prevent sticking
Tissue seared to the forceps tips is ripped away			The forceps can be removed without tissue sticking

Accessories

SILVERGlide Bipolar Forceps

SILVERGlide Bipolar Forceps sterilization trays

- 6700-002-000 Two forceps sterilization tray
- 6700-005-000 Five forceps sterilization tray
- 6700-010-000 10 forceps sterilization tray (forceps not included)

Other SILVERGlide accessories

- 6700-000-002 Polishing cloths
- 6700-000-001 Tip guards
- 6702-000-000 Disposable cable
- 6701-000-000 Reusable cable



Disposable Malis Bipolar Forceps cord and irrigation sets

Rotary cord and tubing sets

- 6790-100-003 Bipolar irrigator integrated tubing and cord set, disposable 6790-100-004 Bipolar irrigator integrated tubing and cord set, disposable, unitized plug

Gravity style cord and tubing sets

- 6790-100-001 Integrated tubing and bipolar cord set, disposable
- 6790-100-002 Integrated tubing and bipolar cord set, disposable, unitized plug

Cord only sets

- 6790-100-005 Bipolar cord set, disposable
- 6790-100-006 Bipolar cord set, disposable, unitized plug



Malis Bipolar Generator and Irrigator system

Malis Bipolar Generator and Irrigator system

- 6700-313-500 Malis Bipolar Generator (includes power cord, does not include foot pedal)
- 6700-313-600 Malis Bipolar Irrigator (does not include power cord)
- 0996-851-011 Universal power cord
- 6700-313-004 Bipolar generator/irrigator interconnect cable
- 6700-313-700 Bipolar generator cart



Irrigator
Optional integrated irrigator

Oversized color screen
for easy comprehension

Programmable user profiles
store information for up to eight surgeons

CUT and COAG knobs facilitate quick setting adjustments

Independent bipolar outputs enhance functionality

DUAL WAVE technology provides confidence during cutting and coagulating

Foot pedals

- 6700-313-003 Single foot pedal (COAG only, does not include required adapter cable)
- 6700-313-304 Single foot pedal adapter cable
- 6700-313-002 Dual foot pedal



Whether you choose reusable or disposable, our non-stick bipolar forceps act as a specialized extension of your hands and eyes.

For more information, contact your Neurosurgical sales representative or call 800 253 3210. You may also learn more at neurosurgical.stryker.com

Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Malis, SILVERglide and Stryker. All other trademarks are trademarks of their respective owners or holders.

NS-GSNPS-BROC-2566610_REV-0
G55/PS
Copyright © 2025 Stryker
Printed in USA

Stryker 1941 Stryker Way
Portage, MI 49002 USA t: 269
323 7700 f: 269 323 2742 toll
free: 800 253 3210

neurosurgical.stryker.com