

CORE Sumex Drill

Your Specialty Is Our Priority.



**Combining the best of our power
and ergonomic designs to give you
a new, high performance drill.**

Sumex Drill

stryker®

Neuro Spine ENT



Zyphr Burs

- Free cutting; efficient, effective bone removal with less effort
- Smoother cuts and reduced chatter
- New flute technology maximizes bone chip clearance and reduces clogging
- New neck design for improved line of sight

Sumex Drill

5400-130-000

- 75000 RPM
- Electronic Torque Feedback (Increases RPM's when the drill is engaged in more demanding tasks)
- Tapered front end (Allows for better line of sight)
- Takes SD/PD Attachments and Maestro Attachments
- Fixed and Steering Duraguards
- Zyphr, Elite and MIS burs
- Complete OR Standardization



Ordering Information

Attachments:

- 5100-120-450 Saber Straight Attachment
- 5100-120-452 Saber Angled Attachment
- 5100-120-470 Saber Long Straight Attachment
- 5100-120-472 Saber Long Angled Attachment
- 5100-120-480 Saber Extra Long Straight Attachment
- 5100-120-482 Saber Extra Long Angled Attachment
- 5100-120-490 Saber MIS Attachment
- 5100-120-492 Saber Angled MIS Attachment
- 5400-210-020 Maestro Straight O Attachment
- 5400-210-050 Maestro Straight M Attachment
- 5400-210-070 Maestro Straight AM Attachment
- 5400-210-080 Maestro R Attachment

- 5400-210-090 Maestro Rx Attachment
- 5400-210-060 Maestro Perforator Chuck
- 5400-210-501 Maestro Adjustable Depth Drill Guide
- 5400-210-011 Maestro Right Angle Attachment

MIS Attachments

- 5100-120-920 SD/PD MIS Medium Straight
- 5100-120-950 SD/PD MIS Long Straight
- 5100-120-922 SD/PD MIS Medium Curved
- 5100-120-952 SD/PD MIS Long Curved

Duraguards:

- 5400-210-057 Maestro Small Fixed Duraguard
- 5400-210-058 Maestro Medium Fixed Duraguard
- 5400-210-059 Maestro Large Fixed Duraguard

- 5400-210-257 Maestro Small Steering Duraguard
- 5400-210-258 Maestro Medium Steering Duraguard
- 5400-210-259 Maestro Large Steering Duraguard
- 5400-210-056 Maestro D-Attachment

Routers:

- 5400-071-057 Small Duraguard Router
- 5400-071-058 Medium Duraguard Router
- 5400-071-059 Large Duraguard Router

Accessories:

- 5400-050-000 Core Console
- 5100-008-000 TPS Footswitch
- 5400-131-000 Hand Switch

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker. Products referenced with [®] designation are registered trademarks of Stryker.

Literature Number: 9100-001-060 Rev. None UnDe/P.S.

Copyright © 2006 Stryker
Printed in USA

Stryker® Neuro Spine ENT
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 269 353 3026
toll free: 800 253 3210

www.stryker.com