

**stryker**

Canadian edition

Sage

**AirTAP™ LC**

**Lift Compatible Patient Repositioning System**

**Advancing mobility**  
with convenience  
and confidence



# Sage AirTAP LC

## Lift Compatible Patient Repositioning System

### Promoting safety and mobility

We understand the importance of early mobility. To achieve your mobility protocols and comply with guidelines, you need a readily available solution that stays under your patient and fits seamlessly into your work flow. That's why we combined our innovative air assist technology with lift compatibility to promote safety and mobility.



**37" x 65"**  
Weight limit: 600 lb / 272 kg

### M<sup>2</sup> Microclimate Body Pad

- Effectively absorbs and locks in moisture to protect patient's skin while allowing air to flow through

### Hi-tech fabric

- Promotes a healthy microclimate for the skin and is compatible with low air loss technology
- Easy to clean surface



### Point of care power switch

- Integrated into workflow to improve efficiency by allowing caregiver to focus on patient at the bedside

### HEPA equipped filter

### Prevalon Air Pump

- Safely fits under most ICU beds
- Hanging clip allows booster pump to travel easily with bed

### Quick connect valve

- Provides an easy, secure connection and a quick release





**30° Body and Anchor Wedge System**

- Helps sustain a 30° side-lying position and reduces pressure by offloading sacrum
- Minimize manual activity as wedges slide into a turn, reducing the need for a log roll
- Minimizes patient migration down the bed
- Reduces the need for boosting, minimizing shear and friction



**Integrated head support**

- Ensures patient's head maintains neutral position while inflated



**Easy roll handles**

- Encourage proper posture for healthcare worker

**Black positioning handles**

- Encourage proper ergonomic position for the healthcare worker when preparing to reposition patient



**Intuitive strap system**

Supine	Reclined	Upright 2	Upright 1

- Enables multiple lifting positions

# 3-in-1 mobility

1

## In-bed positioning

Addresses risk factors of pressure injury

2

## Lateral transfer

Utilizes air-assisted technology to safely transfer patients

3

## Vertical transfer

Enables patients to safely transfer to the bedside chair




# Simple interventions. Extraordinary outcomes.

We are your partner for proven prevention. Our market-leading products solve real healthcare problems and are backed by clinical evidence. Our products allow you to deliver essential patient care with confidence by addressing risk factors that can lead to infections, skin injuries, and caregiver injuries.

**We are driven to solve real problems and make healthcare better for you and your patients.**

## Product **ordering**



### **Sage AirTAP LC Lift Compatible Patient Repositioning System**

- (1) Lift Compatible Glide Sheet
- (1) 30° Body and Anchor Wedge System
- (1) M<sup>2</sup> Microclimate Body Pad

**5 systems/case**  
**Reorder #7287**

### **Sage AirTAP LC Lift Compatible Patient Repositioning System**

#### **Glide Sheet**

- (1) Lift Compatible Glide Sheet
- (1) M<sup>2</sup> Microclimate Body Pad

**10 systems/case**  
**Reorder #7289**

### **30° Body and Anchor Wedge System**

**5 pair/case**  
**Reorder #7295**

### **M<sup>2</sup> Microclimate Body Pad**

36" x 51"  
**30 pads/case (6 bags of 5)**  
**Reorder #7250**



**Prevalon Air Pump  
120V**  
1/case  
**Reorder #7455**



**Prevalon Air Pump  
Cart**  
1 cart/case  
**Reorder #7475**



**Air Pump  
Bed Adapter**  
1/case  
**Reorder #7485**



**HEPA Equipped  
Replacement Filter**  
4/case  
**Reorder #7465**