

Sage

# PrimoFit™

External Urine Management for the Male Anatomy

## A breakthrough in male urine management

PrimoFit provides an easily customizable fit that stays in place and uses continuous suction to effectively manage urinary incontinence while keeping patients comfortable. This first-of-its-kind system serves as an alternative to external collection devices or, when appropriate, indwelling catheters.

PrimoFit's ultra-soft fabric diverts urine away from the skin. Urine is then diverted through the chamber and suctioned into a canister, helping keep the skin dry and addressing a risk factor of incontinence-associated dermatitis (IAD).

**No catheter. No CAUTI.™**

Approximately

**75%**

of hospital-acquired UTIs are associated with an indwelling urinary catheter.<sup>1</sup>

Date: \_\_\_\_\_  
Time: \_\_\_\_\_  
Staff Initials: \_\_\_\_\_  
R6174

### Stays in place

Adhesive pad and base adhesive keeps device in place and allows for easy application and removal



### Customizable fit


Base adhesive eliminates need to use sizing charts and promotes process efficiency

### Helps keep skin dry

Ultra-soft fabric diverts urine away from the skin

### Easy connection

Tiered suction tubing adapter allows for easy connection to hospital suction tubing



Can be left in place for up to **24 hours**



(1) PrimoFit, adult size/package

**40 devices/case**  
**Reorder #5404**

References: 1. Centers for Disease Control and Prevention, Catheter-associated Urinary Tract Infections (CAUTIs). Available at: [https://www.cdc.gov/hai/ca\\_uti/uti.html](https://www.cdc.gov/hai/ca_uti/uti.html). Accessed August 24, 2017.