

# Sage Oral Swab

with Corinz™ Oral Rinse

**REF 6583-X**



## Features and benefits

### Applicator Swab

Applicator Swab

- Helps remove debris and secretions between brushings

Soft foam heads

- Gentle to fragile tissue

## Specifications

Product name	Sage Oral Swab with Corinz™ Oral Rinse
Product code	6583-X
Description	Aids in the removal of debris and oral secretions
Manufacturer accreditation	ISO 13485:2016
Expiration date	24 months from date of manufacture
Length of Suction Swab	15.9 cm
Amount of solution in Corinz Oral Rinse	7 mL
Case quantity	100

### Sodium bicarbonate ingredients

Water (aqua)  
Sodium bicarbonate  
Cellulose gum  
Sodium lauryl sulfate  
Flavour (aroma)  
Sodium saccharin  
Sodium benzoate

### Corinz Oral Rinse ingredients

Water (aqua)  
Glycerin  
Xylitol  
Flavour (aroma)  
Hydroxyethylcellulose  
Potassium sorbate  
Cetylpyridinium chloride  
Sodium saccharin  
Citric acid  
Polysorbate 20  
Polysorbate 80  
Menthol

## Packaging Information

### Contents:

(2) Applicator Swabs with sodium bicarbonate, (1) 7mL Corinz Oral Rinse.

### Indications of Use:

Aids in the removal of debris and oral secretions.

### Instructions for Use:

Burst solution packet with thumbs. Peel Package open. Clean mouth for 1 minute. Instruct to expectorate. Discard according to hospital protocol and local regulations. Repeat as needed/directed.

**Corinz Oral Rinse:** Children under 3 years of age: Consult a dentist or doctor. Children under 12 years of age: Supervise use. Use under the guidance of a healthcare professional. Keep out of reach of children.

### Caution:

Do not allow patient to bite down on oral care tool. Use a bite block if patient has altered levels of consciousness or cannot comprehend commands. Use caution with children and unresponsive individuals. Failure to follow these safety precautions may damage the device and present a choking/aspiration hazard. May not function as intended/potential risk of cross-contamination if device is reused.

Not made with natural rubber latex.

**Applicator Swab:** MDD 93/42/EEC

**Corinz Oral Rinse:** EC Cosmetic Regulation 1223/2009

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local Stryker sales representative.

### Stryker Australia Pty Ltd

8 Herbert Street  
St Leonards NSW 2065  
61 2 9467 1000  
www.stryker.com.au

### Stryker New Zealand Limited

PO Box 17-136, Greenlane  
Auckland 1546, New Zealand  
64 9 573 1890  
www.stryker.com