


## InJector® II Capsule Restoration System

Solving the Limited Working Space Challenge in Hip Arthroscopy



*"The InJector II Capsular Restoration System has been instrumental in my advancement for work in the peripheral compartment. With a T-capsulotomy, I am able to obtain a 180 degree view of the peripheral compartment and address the cam lesion as thoroughly as an open hip dislocation. The InJector II Capsular Restoration System allows me to efficiently close the hip capsule upon completion of the cam osteochondroplasty."*

Shane Nho, M.D., M.S.; Midwest Orthopaedics at Rush, Chicago, IL.

### Product Highlights


- Capsule retraction capability designed to offer complete visualization to the peripheral hip compartment allowing for full access to the cam pathology
- Capsule closure capability is designed to achieve normal capsule anatomy
- Single portal suture passing and retrieval streamlines procedure steps
- Tailhook suture management system allows for highly reproducible stitching even in tough capsule
- Bite depth of 3.5mm (per side) allows for repeatable stitch depth for controlled tissue apposition

## InJector II Ordering Information


Description	Catalog No.
InJector II Handle	CAT01856
InJector II Cartridge	CAT01857
ZipLine #2 Suture 12/box Green/White braid	CAT01585
ZipLine #2 Suture 12/box Black/White braid	CAT01586

For more information about the InJector II Capsule Restoration System, or to request a demo, please call Pivot Medical at 1-866-596-2022

## Anatomic Capsule Restoration



*InJector II enables complete capsule retraction for optimal viewing and resection of the cam pathology.*



*InJector II anatomically restores normal capsule anatomy.*

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: InJector, Stryker. All other trademarks are trademarks of their respective owners or holders.

**Pivot Medical, Inc.**  
A Stryker Sports Medicine Company  
247 Humboldt Court  
Sunnyvale, CA 94089 USA