

# Stronger Than Ever

## Get to the Point.

PEEK Zip Anchor



# PEEK Zip Anchor

## Stronger Than Ever.

---

### **PEEK Zip Anchor**

The PEEK Zip Anchor is fully cannulated and features a trocar-tip inserter to distribute insertion torque along the length of the implant. No tap is needed. The cannulation of the anchor provides a direct channel from the bone marrow to the rotator cuff.

The PEEK Zip Anchor is designed with independent suture eyelets, which are designed to separate each strand of the Force Fiber working suture for easier suture management.

Available in 5.5mm and 6.5mm sizes, with and without half-circle taper needles, the PEEK Zip Anchor is a chosen anchor for rotator cuff and other soft-tissue repairs.

**5.5mm PEEK Zip Anchor**  
3910-200-035

**5.5mm PEEK Zip Anchor with Needles**  
3910-200-030



### 6.5mm PEEK Zip Anchor

3910-200-075



### 6.5mm PEEK Zip Anchor with Needles

3910-200-070



### Features:

---

Trocar-tip inserter to withstand greater torque and for repositioning if necessary.

---

Double-helix thread design for excellent purchase.

---

Separate suture eyelets for better suture management.

---

Available in 5.5mm and 6.5mm sizes.

---

With and without needles.



---

**Joint Replacements**

---

**Trauma, Extremities & Deformities**

---

**Craniomaxillofacial**

---

**Spine**

---

**Biologics**

---

**Surgical Products**

---

**Neuro & ENT**

---

**Interventional Spine**

---

**Navigation**

---

**Endoscopy**

---

**Communications**

---

**Imaging**

---

**Patient Care & Handling Equipment**

---

**EMS Equipment**

325 Corporate Drive  
Mahwah, NJ 07430  
t: 201 831 5000

[www.stryker.com](http://www.stryker.com)

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: **1000-900-580A Rev. 2**  
MS/GS 5/10

Copyright © 2010 Stryker  
Printed in USA