

Suture and tape portfolio

XBraid S, XBraid TT
and Force Fiber




Ordering catalog

A complete suture portfolio

to address your orthopedic needs


Comprehensive portfolio

XBraid and Force Fiber creates a comprehensive suture and tape portfolio with a variety of configurations including specialty options. This portfolio is available in 14 color patterns, various width options, multiple lengths and offered both with and without needles. The XBraid family consists of distinct colors with helical pattern designs that differentiate the braid for maximum visibility.



Strength

XBraid is stronger than the leading competitor with 35 lbf as demonstrated in the knot pull competitive testing below:



Handling characteristics

XBraid's optimal handling characteristics make it smoother during tissue passage and more resistant to cut through.



Low profile knots

XBraid TT (1.4mm) has 228% more surface area and a 59% reduction in knot-stack profile as compared to Arthrex #2 FiberWire.



Arthrex #2 FiberWire Stryker XBraid TT 1.4mm


Cobra Black needles

Cobra Black needles are made with a 300 series stainless steel, giving them at least twice the ductility and superior sharpness than the leading competitor. As compared to stainless steel needles, Cobra Black needles are designed to enhance visibility and reduce glare, especially in blood and smaller surgery sites.



Cobra Black Weave

Cobra Black Weave incorporates helical color patterns designed to enhance visibility.



*Relevant claims and testing data on file at Stryker¹

XBraid independent sutures and tapes

XBraid S, XBraid TT and Force Fiber

XBraid S

Product number	Description
3910900039	#4-0 XBraid S, White/Blue, Suture, BB (Taper, 3/8c, 17mm), 18" Suture
3910900038	#3-0 XBraid S, White/Blue, Suture, BB (Taper, 3/8c, 17mm), 18" Suture
3910900036	#2-0 XBraid S, White/Blue, Suture, CT-3 (Taper, 1/2c, 22mm), 18" Suture
3910900037	#2-0 XBraid S, White/Black, Suture, Non-Needled, 18" Suture
3910900034	#0 XBraid S, White/Blue, Suture, CT-3 (Taper, 1/2c, 22mm), 39" Suture
3910900035	#0 XBraid S, White/Black, Suture, Non-Needled, 39" Suture
3910900033	#1 XBraid S, White/Black, Suture, CT-3 (Taper, 1/2c, 22mm), 46" Suture
3910900023	#2 XBraid S, Blue/Black, Suture, MO-6 (Taper, 1/2c, 26mm), 39" Suture
3910900025	#2 XBraid S, White/Blue, Suture, MO-6 (Taper, 1/2c, 26mm), 39" Suture
3910900026	#2 XBraid S, White/Black, Suture, CT (Taper, 1/2c, 40mm), 39" Suture
3910900027	#2 XBraid S, White/Blue, Suture, Non-Needled, 39" Suture
3910900028	#2 XBraid S, White/Black, Suture, Non-Needled, 39" Suture
3910900029	#2 XBraid S, White/Blue, Suture, ST-1 (Straight Taper Needle, 60mm) on both ends, 39" Suture
3910900030	#2 XBraid S, White/Blue, Suture, CT (Taper, 1/2c, 40mm), 39" Suture
3910900052	#5 XBraid S, White/Black, Suture, CTX (Taper, 1/2c, 48mm), 39" Suture

Force Fiber


Product number	Description
3910900022	#2 Force Fiber, White/Blue, Suture, KC-7 (K Point, 1/2c, 38.9mm), 38" Suture
3910900024	#2 Force Fiber, White/Blue, Suture, C-7 (Taper, 1/2c, 38.9mm), 38" Suture
3910900021	#2 Force Fiber, White/Blue, Suture, Non-Needled, 38" Suture
3910900020	#2 Force Fiber, White/Blue, Suture, HC-5 (Taper, 1/2c, 25.9mm), 38" Suture
CAT01585	#2 Force Fiber, White/Green, Suture, Non-Needled, 36" Suture
CAT01586	#2 Force Fiber, White/Black, Suture, Non-Needled, 36" Suture
3910900051	#5 Force Fiber, White/Blue, Suture, Non-Needled, 38" Suture
3910900050	#5 Force Fiber, White/Blue, Suture, K-60 (3/8c, Sternotomy, 49.2mm), 38" Suture

XBraid TT

Product number	Description
3910900010	1.0mm XBraid TT, Blue/Black, Suture Tape, CT-3 (Taper, 1/2c, 22mm), 28" Suture
3910900017	1.2mm XBraid TT, Violet Co-Braid, 100% UHMWPE
3910900004	1.4mm XBraid TT, White/Blue, Suture Tape, Non-Needled, 46" Suture
3910900013	1.4mm XBraid TT, Black/White, Suture Tape, Non-Needled, 40" Suture
3910900012	1.4mm XBraid TT, White/Blue, Suture Tape, MO-6 (Taper, 1/2c, 26mm), 40" Suture
3910900011	1.4mm XBraid TT, Blue/Black, Suture Tape, CT (Taper, 1/2c, 40mm), 40" Suture
3910900002	1.4mm, White/Blue, Suture Tape, Non-Needled, 40" Overall Length, 9" tape working length
3910900003	1.4mm, White/Black, Suture Tape, Non-Needled, 40" Overall Length, 9" tape working length
3910900006	1.4mm XBraid TT, Blue/Black, Suture Tape, MO-6 (Taper, 1/2c, 26mm), 40" Suture
3910900005	1.4mm XBraid TT, White/Blue, Suture Tape, CT (Taper, 1/2c, 40mm), 40" Suture
3910900009	1.4mm XBraid TT, White/Blue, Suture Tape, ST-1 (Straight Taper Needle, 60mm), on both ends, 40" Suture
3910900018	2.0mm XBraid TT, Violet Co-Braid, 100% UHMWPE
3910900019	2.0mm XBraid TT, Violet, 100% UHMWPE
3910900015	2.2mm XBraid TT, White/Blue, Suture Tape, MO-6 (Taper, 1/2c, 26mm), 40" Suture
3910900016	2.2mm XBraid TT, White/Black, Suture Tape, CTX (Taper, 1/2c, 48mm), 40" Suture
3910900007	2.2mm XBraid TT, White/Blue, Suture Tape, CTX (Taper, 1/2c, 48mm), 40" Suture
3910900008	2.2mm XBraid TT, White/Black, Suture Tape, MO-6 (Taper, 1/2c, 26mm), 40" Suture



XBraid specialty sutures and tapes



Passing Loop


XBraid Passing Loops allow for multiple sutures to be shuttled with one pass and facilitate a luggage tag technique.

Product number	Description
3910900042	#2-0 XBraid S, White/Blue, Passing Loop Suture, Non-Needled, 24" Suture w/ 1.5" Loop
3910900041	#0 XBraid S, White/Blue, Passing Loop Suture, Non-Needled, 24" Suture w/ 1.5" Loop
3910900040	#2 XBraid S, White/Blue, Passing Loop Suture, Non-Needled, 24" Suture w/ 1.5" Loop
3910900044	1.0mm XBraid TT, Blue/Black, Passing Loop Suture Tape, Non-Needled, 24" Suture w/ 7" Loop
3910900043	1.4mm XBraid TT, White/Blue, Passing Loop Suture Tape, Non-Needled, 24" Suture w/ 5" Loop

XPass

XBraid XPass is a 50" suture that provides a 12" stiffened end that is designed to pass through cannulation or needles.


Product number	Description
3910900072	#2-0 XBraid S, White/Blue, XPass Suture, Non-Needled, 50" Suture (one end stiffened 12")
3910900070	#2 XBraid S, White/Blue, XPass Suture, Non-Needled, 50" Suture (one end stiffened 12")
3910900071	#2 XBraid S, White/Black, XPass Suture, Non-Needled, 50" Suture (one end stiffened 12")



XPass with Loop

XBraid XPass with Loop provides a looped suture with a stiffened tip to facilitate passage through bone tunnels

Product number	Description
3910900076	#2 XBraid S, White/Blue, XPass Suture w/ Loop, Non-Needled, 26" Suture (one end stiffened 12", other end .75" Loop)



Whip Stitch


XBraid Whip Stitch has tape and suture options, ideal for multi-strand tendon repair. Loops are available in various sizes and colors including a smaller 10" option up to a 20" option.

Product number	Description
3910900063	#4-0 XBraid S, White/Blue, Whip Stitch Suture, BB (Taper, 3/8c, 17mm), 10" Loop length
3910900062	#2-0 XBraid S, White/Blue, Whip Stitch Suture, CT-3 (Taper, 1/2c, 22mm), 20" Loop length
3910900061	#0 XBraid S, White/Blue, Whip Stitch Suture, CT-3 (Taper, 1/2c, 22mm), 20" Loop length
3910900060	#2 XBraid S, White/Blue, Whip Stitch Suture, ST-1 (Straight Taper Needle, 60mm) 20" Suture w/ 18" Loop
3910900064	1.4mm XBraid TT, White/Blue, Whip Stitch Suture Tape, ST-1 (Straight Taper Needle, 60mm), 20" Suture w/ 18" Loop











SharpShooter

XBraid SharpShooter includes 13" double armed long trocar meniscal needles to provide an excellent option for standard inside-out meniscal repair procedures.


Product number	Description
HS4702	#2-0 SharpShooter Suture, White/Blue, 13" Straight Trocar Needle on both ends, 30" Suture
HS4703	1.0mm SharpShooter Suture Tape, Blue/Black, 13" Straight Trocar Needle on both ends, 30" Suture w/ 5" Suture Tape
HS4704	#2-0 SharpShooter Suture, White UHMWPE Suture, 13" Straight Trocar Needle on both ends, 36.5" Suture
PB4704	#2-0 SharpShooter Suture, Green Polyester Suture, 13" Straight Trocar Needle on both ends, 36" Suture



Needle options

Product description	Needle image
BB, 3/8C, 17mm	
CT-3, 1/2C, 22mm	
HC-5, 1/2C, 25.9mm	
MO-6, 1/2C, 26mm	
C-7, 1/2C, 38.9mm	
KC-7, 1/2 C, 38.9mm	
CT, 1/2C, 40mm	
CTX, 1/2 C, 48mm	
K-60, 3/8C, 49.2mm	
ST-1, Straight, 60mm	

Sample color options



- | | | |
|---|--|--|
| A = 1.2mm XBraid TT
Violet Co-Braid | F = 1.5mm XBraid TT
White Blue | K = 1.4mm XBraid TT
Blue/Black* |
| B = 2.0mm XBraid TT
Violet | G = 1.4mm XBraid TT
White/Blue | L = 1.4mm XBraid TT
White/Black* |
| C = 2.2mm XBraid TT
White/Blue | H = #2 Force Fiber
White/Green | M = #2 XBraid S
White/Black* |
| D = #5 Force Fiber
White/Blue | I = #2 XBraid S,
White/Blue | N = #2 XBraid S
Blue/Black* |
| E = #2 Force Fiber
White/Blue | J = 1.4mm XBraid TT
White/Black* | *=Cobra Black Weave |

References

- Riverpoint Medical Publication PUB375-01, 2020 UHMWPE Informational Datasheet.

Sports Medicine

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/ or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets.

Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: XBraid, Force Fiber and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

Cobra Black is a registered trademark of RiverPoint Medical.

All features of competitive products were obtained from publicly available information.

1000903307 Rev B

Copyright © 2020 Stryker

Stryker
Sports Medicine
5670 Greenwood Plaza Blvd.
Ste. 200
Greenwood Village, CO 80111

t: 866 596 2022
www.sportsmedicine.stryker.com