

Neptune® SafeAir® **Smoke Evacuation Portfolio**



Ordering information



Pencil – push button switch	
Part number	Description
0703-046-000	Integrated smoke evacuation pencil, coated 70mm blade, push button switch
0703-046-001	Integrated smoke evacuation pencil, non-coated 70mm blade, push button switch
0703-046-002	Non-sterile integrated smoke evacuation pencil, coated 70mm blade, push button switch
0703-046-003	Non-sterile integrated smoke evacuation pencil, non-coated 70mm blade, push button switch

Pencil – rocker switch	
0703-047-000	Integrated smoke evacuation pencil, coated 70mm blade, rocker switch
0703-047-001	Integrated smoke evacuation pencil, non-coated 70mm blade, rocker switch
0703-047-002	Non-sterile integrated smoke evacuation pencil, coated 70mm blade, rocker switch
0703-047-003	Non-sterile integrated smoke evacuation pencil, non-coated 70mm blade, rocker switch

Telescopic pencil – push button switch	
0703-046-004	Integrated telescopic smoke evacuation pencil; coated 70mm blade; push button switch
0703-046-005	Integrated telescopic smoke evacuation pencil; uncoated 70mm blade; push button switch
0703-046-006	Integrated telescopic smoke evacuation pencil; coated 70mm blade; push button switch, non-sterile*

Telescopic pencil – rocker switch	
0703-047-004	Integrated telescopic smoke evacuation pencil; coated 70mm blade; rocker switch
0703-047-005	Integrated telescopic smoke evacuation pencil; uncoated 70mm blade; rocker switch
0703-047-006	Integrated telescopic smoke evacuation pencil; coated 70mm blade; rocker switch, non-sterile*

*Non-sterile options for procedure packs available through distributors

Suction sleeve		
0703-005-065	Suction sleeve	65mm
0703-005-070	Suction sleeve	70mm
0703-005-125	Suction sleeve	125mm
0703-005-165	Suction sleeve	165mm

SafeAir Compact smoke evacuator	
0703-045-000	SafeAir Compact smoke evacuator, stand alone
0703-040-000	ULPA filter, compact smoke evacuator
0703-045-100	Smoke tubing for canister/fluid trap
0703-045-001	Foot pedal, compact smoke evacuator
0703-040-001	Fluid trap

Needle		
Part number	Description	Length
0703-007-070	Needle electrode	70mm
0703-007-125	Needle electrode	125mm

Blade		
0703-070-000	Blade electrode coated	70mm
0703-070-001	Blade electrode	70mm
0703-070-002	Blade electrode coated and insulated, 5mm exposed	70mm
0703-070-003	Blade electrode coated and insulated, 2mm exposed	70mm
0703-070-005	Blade electrode coated square	70mm
0703-125-000	Blade electrode coated	125mm
0703-125-001	Blade electrode	125mm
0703-125-002	Blade electrode coated and insulated	125mm
0703-165-000	Blade electrode coated	165mm
0703-165-001	Blade electrode	165mm
0703-165-002	Blade electrode coated and insulated	165mm

Ball		
0703-120-003	Ball electrode 3mm diameter	120mm
0703-120-005	Ball electrode 5mm diameter	120mm

Loop and conization				
Part number	Description	Width	Depth	Length
0703-213-015	Conization electrode	13mm	15mm	120mm
0703-213-020	Conization electrode	13mm	20mm	120mm
0703-216-008	Conization electrode	16mm	8mm	120mm
0703-216-018	Conization electrode	16mm	18mm	120mm
0703-220-020	Conization electrode	20mm	20mm	120mm
0703-310-010	Loop t-bar electrode	10mm	10mm	120mm
0703-315-012	Loop t-bar electrode	15mm	12mm	120mm
0703-320-015	Loop t-bar electrode	20mm	15mm	120mm
0703-320-020	Loop t-bar electrode	20mm	20mm	120mm
0703-620-015	Loop t-bar electrode	20mm	15mm	60mm
0703-620-020	Loop t-bar electrode	20mm	20mm	60mm
0703-720-020	Loop u-bar electrode	20mm	20mm	120mm

The SafeAir Smoke Evacuation Pencils are also compatible with Stryker’s full offering of the Colorado MicroDissection needles.

2-in-1 solution

The SafeAir Smoke Evacuation Pencils provide a 2-in-1 solution that feels natural to use. Utilizing monopolar energy, they enable smoke evacuation, cutting and coagulation via a range of electrodes.

The following components come standard with your initial purchase, outfitting you with the basics needed to perform a spectrum of procedures.

- The SafeAir Smoke Evacuation Pencil with your choice of push button or rocker switch
- 70mm blade electrode, coated or uncoated based on your preference
- Suction sleeve
- Protective holster**
- Tubing and drape clip**

**Not available with the SafeAir Telescopic Smoke Evacuation Pencil

Try it today!

To learn more about the Neptune SafeAir Smoke Evacuation Portfolio, to request a product demo or learn about Stryker's Smoke Free O.R. Program, please contact your Surgical sales representative or call **800 253 3210**. To learn more about the hazards of surgical smoke and our commitment to smoke evacuation, visit neptunewastemanagement.com/SafeAir.

Surgical Technologies

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Colorado, Neptune, SafeAir and Stryker. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.