

SurgiCount® Safety-Sponge® System



Product
catalog



SurgiCount

Safety-Sponge System

SurgiCount tablet

Since 2009, **more than 550** U.S. healthcare facilities have trusted the SurgiCount Safety-Sponge System, using **200 million** of our sponges in **11 million** procedures with **zero** retained surgical sponges. Here's a look at the time-tested design and capabilities that lead to this success – plus the next-generation features that elevate your benefits to a new level.

- Technology that grows with you
- Expandable technology platform applicable well into the future
- Secure cloud-based server/storage
- Intuitive, user-friendly tablet interface
- WiFi capability with real-time data transfer to server (or send via hard-wired docking)
- Admin portal for password protected, VPN access from any onsite or remote computer
- Workflow simplicity
- Tablet self-reconciliation can be captured through SurgiCount360 software and customized for post-op compliance/reporting
- Adds only one second to mean count time (mean, 5 seconds to 4 for a non-barcoded sponge), but no increase in overall case length since activities are done in parallel
- One-step set up and turnover
- ID-based count reconciliation can eliminate added steps for detection scanning or wandling
- Doesn't require added capital or disposable costs for detection mats, detection wands and sterile wand sleeves; and avoids associated medical waste





SurgiCount

Safety-Sponge System

SurgiCounter kit (0694-001-000) includes



SurgiCount Tablet kit
0694-001-010



Tablet pole mount
0694-001-050



Tablet dock
0694-001-040



Tablet battery (x3)
0694-001-030



Tablet power supply
0694-001-080

Other tablet accessories



Tablet dock wall mount
0694-001-060



Tablet wall mount
0694-001-070



SurgiCount360 Software Application

With expanded functionality, SurgiCount360 Software provides accurate post-operative fully compliant documentation—a vital part of the Safety-Sponge System. By downloading information captured intra-operatively, the SurgiCount360 Software allows access to both individual procedure reports and aggregate data to provide auditable evidenced based outcomes and enhanced patient safety.

- Options allow viewing over a customized time frame and a variety of purposes including
 - Post-operative evidence that all individual SurgiCount Safety-Sponge products used have been successfully accounted for outside of the patient
 - Monitoring adherence to the facility's sponge/lap/towel count procedure
 - Providing detailed SurgiCount Safety-Sponge product usage data
 - Designed for both internal use and third party auditors
- Easy-to-use interface provides simple data access and review
- Email alert function provides information in situations with potential compliance or patient safety issues



Gauze

SM-2212-6S

Gauze

- 2x2
- 12 ply
- Qty per pack - 10
- Packs per case - 128
- Master tagged
- Xray



SM-3312-6S

Gauze

- 3x3
- 12 ply
- Qty per pack - 10
- Packs per case - 128
- Master tagged
- Xray



SM-3312-FS

Gauze

- 3x3
- 12 ply
- Qty per pack - 10
- Packs per case - 128
- Master tagged
- Xray



Gauze

SM-4416-7S*

Gauze

- 4x4
- 16 ply
- Qty per pack - 10
- Packs per case - 128
- Master tagged
- Double xray



SM-4416-9S

Gauze

- 4x4
- 16 ply
- Qty per pack - 10
- Packs per case - 128
- Master tagged
- Double wavy xray



SM-4432-7S

Gauze

- 4x4
- 32 ply
- Qty per pack - 10
- Packs per case - 128
- Master tagged
- Double xray



SM-8412-7S*

Gauze

- 8x4
- 12 ply
- Qty per pack - 10
- Packs per case - 48
- Master tagged
- Double xray



Gauze

SM-8416-7S

Gauze

- 8x4
- 16 ply
- Qty per pack - 10
- Packs per case - 48
- Master tagged
- Double xray



SM-8424-7S

Gauze

- 8x4
- 24 ply
- Qty per pack - 10
- Packs per case - 48
- Master tagged
- Double xray



Lap sponges

SX-0418-PS*

SAFE-T lap sponge

- 4x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Enhanced xray visibility



SM-0418-PS*

Lap sponge

- 4x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray



SR-0418-PS

Lap sponge

- 4x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Enhanced xray
- Plastic ring



SX-0836-PS*

SAFE-T lap sponge

- 8x36
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Enhanced xray visibility



Lap sponges

SM-0836-PS

Lap sponge

- 8x36
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray



SR-0836-PS

Lap sponge

- 8x36
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray
- Plastic ring



SX-1212-PS*

SAFE-T lap sponge

- 12x12
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Enhanced xray visibility



SM-1212-PS*

Lap sponge

- 12x12
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray



Lap sponges

SR-1212-PS

Lap sponge


- 12x12
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray
- Plastic ring between layers



SM-1726-PS

Lap sponge

- 17x26
- Qty per pack - 5
- Packs per case - 20
- Master tagged
- Xray



SX-1818-PS*

SAFE-T lap sponge

- 18x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Enhanced xray visibility



SM-1818-PS*

Lap sponge

- 18x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray



Lap sponges

SR-1818-PS

Lap sponge

- 18x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray
- Plastic ring between layers



SM-1818-NS

Lap sponge

- 18x18
- Qty per pack - 5
- Packs per case - 40
- Master tagged
- Xray
- No loop



SM-1836-PS

Lap sponge

- 18x36
- Qty per pack - 3
- Packs per case - 20
- Master tagged
- Xray



Vaginal lap sponges

SD-0418-VS

Vaginal lap sponge

- 4x18
- Qty per pack - 1
- Packs per case - 128
- Xray



SD-0836-VS

Vaginal lap sponge


- 8x36
- Qty per pack - 1
- Packs per case - 128
- Xray



SV-0418-US

Vaginal lap sponge

- 4.5x18
- Qty per pack - 1
- Packs per case - 40
- 16 ply
- Xray
- Not pre-washed



SV-0433-US

Vaginal lap sponge

- 4.5x33
- Qty per pack - 1
- Packs per case - 40
- 8 ply
- Xray
- Not pre-washed



Vaginal lap sponges

SV-0436-PS

Vaginal lap sponge

- 4x36
- Qty per pack - 1
- Packs per case - 40
- 4 ply
- Xray



SV-0836-PS

Vaginal lap sponge

- 8x36
- Qty per pack - 1
- Packs per case - 40
- 4 ply
- Xray



SV-0236-6S

Vaginal gauze roll

- 2x36
- Qty per pack - 1
- Packs per case - 40
- 4 ply
- Xray



SV-0272-6S

Vaginal gauze roll

- 2x72
- Qty per pack - 1
- Packs per case - 40
- 4 ply
- Xray



O.R. towels

SM-6002-BS

O.R. towel

- 17x26
- Qty per pack - 2
- Packs per case - 40
- Blue
- Master tagged
- Xray with loop



SM-6002-GS

O.R. towel


- 17x26
- Qty per pack - 2
- Packs per case - 40
- Green
- Master tagged
- Xray with loop



SM-6002-WS

O.R. towel

- 17x26
- Qty per pack - 2
- Packs per case - 40
- White
- Master tagged
- Xray with loop



SM-6005-BS

O.R. towel

- 17x26
- Qty per pack - 5
- Packs per case - 16
- Blue
- Master tagged
- Xray with loop



O.R. towels

SM-6005-GS

O.R. towel


- 17x26
- Qty per pack - 5
- Packs per case - 16
- Green
- Master tagged
- Xray with loop



SM-6005-WS*

O.R. towel

- 17x26
- Qty per pack - 5
- Packs per case - 16
- White
- Master tagged
- Xray with loop



Ordering information

For questions or to place your order today,
please contact your Stryker Surgical sales
representative or call

800 253 3210

Stryker Instruments

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: SAFE-T Lap, Safety-Sponge, Safety-Sponge System, SurgiCount, SurgiCount360 and Stryker.

All other trademarks are trademarks of their respective owners or holders.

Literature Number: 9100-002-539-Rev E
Copyright © 2018 Stryker

Stryker Instruments
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210
stryker.com