


Tray layout card – T2 Alpha Tibia Indication Tray

T2 Alpha Tibia Nailing System



Base tray



Ref #	Description	Qty.	Ref #	Description	Qty.
A 2355-4000	Tibia Indication Tray		4 2353-4104	Nail Holding Screw Tibia / Femur PF	2
1 2351-0255	External Compression Device Tibia	1	Optional		
2 2351-4107	Targeting Arm Tibia	1	5 2355-5030	Drills and Pins Insert	1
3 2353-4103	Nail Adapter Tibia	1	6 2355-4010	SPI Insert	1

Tray layout card – T2 Alpha Tibia Indication Tray



Opening 11.5mm Insert



Ref #	Description	Qty.	Ref #	Description	Qty.
1 2351-6112	Opening Reamer	1	4 2351-6312	Opening Reamer Trocar	1
2 2351-6000	Opening Reamer Handle	1	Optional		
3 2351-6212	Opening Reamer Sleeve	1	5 2351-0145	Quick-Lock T-Handle, Modified Trinkle	1

Tray layout card – T2 Alpha Tibia Indication Tray



SPI Insert layout and content (optional)



Ref #	Description	Qty.	Ref #	Description	Qty.
1	1806-1412	1	6	1806-1415	1
2	1806-1408	1	7	2353-4113	1
3	1806-1414	1	8	1806-1416	1
4	2353-4114	1	9	1806-1400	1
5	1806-1409	1	10	1806-1417	2

Tray layout card – T2 Alpha Tibia Indication Tray

Drills and pins insert layout and content (optional)



Ref #	Description	Qty.	Ref #	Description	Qty.
	2351-3340S* Threaded Recon K-Wire; 3.2 × 400mm				
1	2351-3240S* Drill-Tip Recon K-Wire; 3.2 × 400mm	4	5	2351-5510 Counterbore Drill Ø5, Long	
	2351-3028S* K-Wire; 3 × 285mm		6	2351-4213S* Freehand Drill; 4.2 × 130mm	2
2	2351-4236S* Locking Drill; 4.2 × 360mm	2		2351-3513S* Freehand Drill; 3.5 × 130mm	
	2351-3536S* Locking Drill; 3.5 × 360mm		7	2351-5515 Counterbore Drill Ø5, Manual	1
3	2351-5500 Counterbore Drill Ø5, Short	1			
4	2351-3518S* Freehand Drill; 3.5 × 185mm	2			
	2351-4218S* Freehand Drill; 4.2 × 185mm				

**sold as prepackaged sterile*

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area. Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, T2 Alpha. All other trademarks are trademarks of their respective owners or holders.

Content ID: T2-TLC-4, 09-2018
Copyright © 2018 Stryker



Manufacturer:
Stryker Trauma GmbH
Prof.-Küntschner-Str. 1-5
24232 Schönkirchen
Germany

stryker.com