


stryker®

Trauma & Extremities

VariAx[®] 2 Wrist Fusion

System Components Guide

Stryker Plating



Contents

	Page
1. VariAx 2 Wrist Fusion Tray	3
Top Level & 2 Level Tray Frame	3
Drawer 1	4

TOP LEVEL & 2 LEVEL TRAY FRAME

CONSISTING OF:



REF	Description
940090	Wrist Fusion Tray (including Labeling Clip)
940093	Labeling Clip Wrist Fusion Tray (Spare Part)

VARIAX 2 INSTRUMENTS



REF	Description
703920	Handle, Large
703921	Handle, Medium
703923	Handle Insert, AO, cannulated
703922	Handle Insert, AO, Ratchet, cannulated
703663	Screwdriver Blade, AO, T8, Self Retaining
703880	Screwdriver Blade, AO, T10, Self Retaining
703927	Joystick for T8 screw holes
703928	Joystick for T10 screw holes
705170	Depth Gauge 2.4mm - 3.5mm VariAx Screws
703684	Drill Guide, 2.0mm Drill, Comp/Polyaxial T8
703882	Drill Guide, 2.6mm Drill, Comp/Polyaxial T10
703884	Drill Guide, For 2.7mm independent lag screw T8
703686	Drill Guide, For 3.5mm independent lag screw T10
703896	Drill Bit, AO, Ø2.0mm x 135mm, Scaled
703691	Drill Bit, AO, Ø2.6mm x 135mm, Scaled
703897	Overdrill, AO, Ø2.7mm x 122mm
703694	Overdrill, AO, Ø3.5mm x 122mm
45-80040	Countersink for T8/T10 screws, AO
703898	Tap, AO, For 3.5mm x 122mm Screws
703899	Tap, AO, For 2.7mm x 102mm Screws
703900	Tap, AO, For 2.4mm x 102mm Screws
703886	SpeedGuide™, for 3.5mm Screws, T10 (L = 30mm)
703888	SpeedGuide™, for 2.4/2.7mm Screws, T8 (L = 30mm)
703891	SpeedGuide™ Drill, AO, Ø2.0mm (L = 30mm)
703894	SpeedGuide™ Drill, AO, Ø2.6mm (L = 30mm)


• All non-sterile drills, taps and countersink may be ordered sterile by placing an "S" at the end of the REF Number



VariAx 2 Wrist Fusion Tray

DRAWER 1

WRIST FUSION PLATES, INSTRUMENTS & SCREWS



REF	Description
1500-0005	Drawer For Modules & Screw Racks
940096	Labeling Clip Wrist Fusion Drawer
940092	Wrist Fusion Implant & Instrument Insert
940094	Wrist Fusion Screw Rack T8
940095	Wrist Fusion Screw Rack T10

WRIST FUSION PLATES



REF	Description
629579	Standard Bend Plate
629589	Short Bend Plate
629590	Straight Plate

WRIST FUSION PLATE TRIALS



REF	Description
705305	Standard Bend Plate Trial
705306	Short Bend Plate Trial
705307	Straight Plate Trial

INSTRUMENTS



REF	Description
703938	Bending Iron for VariAx Plates
705291	Periosteal Elevator 3.0mm round edges, straight blade
705292	Periosteal Elevator 3.0mm round edges, curved blade (optional)
703818	2mm K-Wire with Stop (pack of 5)
390164	K-Wire 1.6 x 150mm (pack of 10)

3.5MM LOCKING SCREW

Titanium REF	Length mm
	657312 12
	657314 14
	657316 16
	657318 18
	657320 20
	657322 22
	657324 24
	657326 26
	657328 28
	657330 30
	657332 32

3.5MM NON LOCKING SCREW

Titanium REF	Length mm
	657412 12
	657414 14
	657416 16
	657418 18
	657420 20
	657422 22
	657424 24
	657426 26
	657428 28
	657430 30
	657432 32

• All non-sterile plates, screws, and k-wires may be ordered sterile by placing an "S" at the end of the REF Number

**2.7MM LOCKING SCREWS**

Titanium REF	Length mm
656310	10
656312	12
656314	14
656316	16
656318	18
656320	20
656322	22
656324	24

**2.7MM NON LOCKING SCREWS**

Titanium REF	Length mm
656410	10
656412	12
656414	14
656416	16
656418	18
656420	20
656422	22
656424	24

**2.4MM LOCKING SCREW**

Titanium REF	Length mm
656010	10
656012	12
656014	14
656016	16
656018	18
656020	20
656022	22
656024	24

**2.4 NON LOCKING SCREWS**

Titanium REF	Length mm
656110	10
656112	12
656114	14
656116	16
656118	18
656120	20
656122	22
656124	24

• All non-sterile screws may be ordered sterile by placing an "S" at the end of the REF Number

Notes

Notes

Reconstructive

Hips
Knees
Trauma & Extremities
Foot & Ankle
Joint Preservation
Orthobiologics & Biosurgery

MedSurg

Power Tools & Surgical Accessories
Computer Assisted Surgery
Endoscopic Surgical Solutions
Integrated Communications
Beds, Stretchers & EMS
Reprocessing & Remanufacturing

Neurotechnology & Spine

Craniomaxillofacial
Interventional Spine
Neurosurgical, Spine & ENT
Neurovascular
Spinal Implants

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for Cleaning and Sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker, VariAx. All other trademarks are trademarks of their respective owners or holders.

The SmartLock Technology is patented by Professor Dietmar Wolter, Hamburg Germany.

The products listed above are CE marked.

Content ID: VAX-SCS-6, 6-2015

Copyright © 2015 Stryker



Manufactured by:

Stryker GmbH
Bohnackerweg 1
2545 Selzach
Switzerland

www.stryker.com