

Insignia™ Hip Stem

Founded on Stryker's signature principles of implant fit, Insignia is your collared, **metaphyseal-filling** hip stem designed specifically for muscle-sparing approaches.



*Optimization subject to particular design constraints

Tri-Stage™ Broach

Extraction



M/L: Extraction teeth facilitate cancellous bone removal for enhanced cortical fit¹

Compaction



A/P: Compaction teeth enable bone preservation and initial stability

Distally Cutting



Distal: Distal diamond cutting teeth help remove diaphyseal bone and are oversized by 0.5mm to promote proximal fit



Ordering information

Insignia™ Hip Stem

	Standard	High Offset
Size 0	7000-5500	7000-6600
Size 1	7000-5501	7000-6601
Size 2	7000-5502	7000-6602
Size 3	7000-5503	7000-6603
Size 4	7000-5504	7000-6604
Size 5	7000-5505	7000-6605
Size 6	7000-5506	7000-6606
Size 7	7000-5507	7000-6607
Size 8	7000-5508	7000-6608
Size 9	7000-5509	7000-6609
Size 10	7000-5510	7000-6610
Size 11	7000-5511	7000-6611

Insignia Hip Stem material composition is Ti-6Al-4V ELI as per specification ASTM F136 (wrought).

References:

1. Insignia Design Verification Memo D0000097336.

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate the breadth of Stryker's product offerings.

A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: [insert names of appropriate trademarks in alphabetical order]. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

INSIGN-SS-1_31783

Copyright © 2021 Stryker