
stryker

VersiQuad

Quadriceps harvesting system



Engineered for sharp cuts, personalized options and clear visibility[1,2]

Clean, customized and clear cuts[1,3]

The VersiQuad System is a sharp harvesting system designed for **increased visualization, customization and a minimally invasive harvest**[1,3].


Featuring both partial and full-thickness parallel blades and eight size options, it allows you to tailor graft size to your preference[1]. The VersiQuad Harvesting System delivers comprehensive, start-to-finish support throughout the quadriceps autograft harvesting process.


Parallel blade


- The calibrations on the handle provide estimation of total graft length[1]
- **Depth options** allow for surgeon choice of full thickness harvesting or partial thickness harvesting technique[1]
- Customization of graft size with **8 size options**
- Single-use sharp blade with reusable handle[3]
- Radiused blade allows for cutting in a push or pull motion[3]
- The disposable blades are designed to be easily inserted using the BS EN 27740:1992 instrumentation standard[1]

Tendon cutter


- 12 x 12 cutting window to accommodate larger grafts and bone blocks[2]
- Sharp, single use disposable blade with reusable handle designed with sharpness of cut for each use[4]
- Locking lever to aid in prevention of premature termination of graft[3]
- Graft length estimations on handle[1]







Visualization for surgical performance



Halo Retractor

- Accommodates a minimally invasive incision[2]
- 10mm incremental calibrations to estimate graft length[5]
- Accommodates both a vertical or horizontal incision[1]
- Reduces the need for general surgical retracting instruments[1]



Invuity Lighted Retractor

- Designed with a reusable light cord and retractor handle to facilitate visualization of the graft harvest site[6]
- Compatible with commercially available Xenon* and LED light sources[7]

Featured products

Product number	Description
234-112-500	Tendon Cutter
234-112-600	Parallel Blade Handle
234-112-012	Tendon Cutter Blade
234-112-108	Parallel Blade 8mm Width x 10mm Thickness
234-112-109	Parallel Blade 9mm Width x 10mm Thickness
234-112-110	Parallel Blade 10mm Width x 10mm Thickness
234-112-607	Parallel Blade 7mm Width x 6mm Thickness
234-112-608	Parallel Blade 8mm Width x 6mm Thickness
234-112-609	Parallel Blade 9mm Width x 6mm Thickness
234-112-610	Parallel Blade 10mm Width x 6mm Thickness
234-112-611	Parallel Blade 11mm Width x 6mm Thickness
234-112-100	Skin Retractor
C40XN	Invuity Lighted Retractor 40mmx12mm Extra Narrow [Heaney]
C3550ANXN	Invuity Army Navy Lighted Retractor Extra Narrow
104008	Narrow Invuity Photon Guide Illuminator
FC1S	Single fiber optic cable, ACMI 10 ft
INV-LLS-2000HPAK	LED Turbo Light Source
234102067	ProCinch Adjustable Loop Implant QuadCinch with button
234102068	ProCinch Adjustable Loop Implant QuadCinch no button
0234100001	ProCinch slotted button concave round 11mm
0234100002	ProCinch slotted button concave round 14mm
0234100003	ProCinch slotted button concave round 20mm
0234100004	ProCinch slotted button flat 8 X 12mm
0234100005	ProCinch slotted button flat 14mm

References

1. TR24134
2. DHD15249
3. TR23709
4. DHD15164 Rev C
5. P57075
6. IFU - 700000772135
7. IFU 700000770163

*Only use Xenon Light Sources that are less than 300 watts at less than 50% output setting, with a minimum distance of 2cm from the surrounding tissue. Refer to the Illuminated Retractor System Instructions for Use

Sports Medicine

This information is intended solely for the use of healthcare professionals. A healthcare professional must always rely on their own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using the product.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker Representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use, or have applied for the following trademarks or service marks: VersiQuad, Stryker. All other trademarks are trademarks of their respective owners or holders.

Stryker
Sports Medicine
5670 Greenwood Plaza Blvd, Ste 200
Greenwood Village, CO 80111

www.stryker.com

SM-KNE-BROC-3655143

Copyright © 2026 Stryker